

# ELECTRIC MODELS ARE ALSO AVAILABLE



## Model 25X (25 qt/24 liter)

**Model 25X-120:** 120 Volt, 50/60 Hz 1050 watts/8.75 amps

**Model 25X-240:** 240 Volt, 50/60 Hz 1050 watts/4.38 amps

<b>Gross Capacity</b>	<b>25 qt / 24 liter</b>
Overall Height	16 <sup>3</sup> / <sub>4</sub> " / 42.5cm
Bottom Height	12 <sup>1</sup> / <sub>4</sub> " / 31.2cm
Inside Diameter	12 <sup>5</sup> / <sub>8</sub> " / 32.1cm
Unit Weight	26 lbs. / 11.8kg.

### Inner Container No. 2156

Height	8 <sup>1</sup> / <sub>2</sub> " / 21.6cm
Diameter	11 <sup>1</sup> / <sub>8</sub> " / 28.3cm
Circumference	35 <sup>7</sup> / <sub>8</sub> " / 91.1cm
Capacity	14.5 qt / 13.7 liter
Volume	835in <sup>3</sup> / 13,688cm <sup>3</sup>

Carton Dimensions	19 <sup>1</sup> / <sub>2</sub> " x 17 <sup>1</sup> / <sub>2</sub> " x 19"
	56.5cm x 44.5cm x 48.3cm

Shipping Weight	30 lb. / 13.6kg.
Unit Pack: 1	Cube: 3.75

<b>Optional No. 2180 Support Base</b>	<b>2" / 5cm high</b>
Outside Diameter	12 <sup>3</sup> / <sub>4</sub> " / 32.4cm
Inside Diameter	12 <sup>3</sup> / <sub>8</sub> " / 31.4cm
Elevates Sterilizer Above Surface	1 <sup>3</sup> / <sub>8</sub> " / 3.5cm



## Model 50X (25 qt/24 liter)

**Model 50X-120:** 120 Volt, 50/60 Hz 1650 watts/13.75 amps

**Model 50X-240:** 240 Volt, 50/60 Hz 1650 watts/6.88 amps

<b>Gross Capacity</b>	<b>25 qt / 24 liter</b>
Overall Height	16 <sup>3</sup> / <sub>4</sub> " / 42.5cm
Bottom Height	12 <sup>1</sup> / <sub>4</sub> " / 31.2cm
Inside Diameter	12 <sup>5</sup> / <sub>8</sub> " / 32.1cm
Unit Weight	29 lbs. / 13.2kg.


### Inner Container No. 2156

Height	8 <sup>1</sup> / <sub>2</sub> " / 21.6cm
Diameter	11 <sup>1</sup> / <sub>8</sub> " / 28.3cm
Circumference	35 <sup>7</sup> / <sub>8</sub> " / 91.1cm
Capacity	14.5 qt / 13.7 liter
Volume	835in <sup>3</sup> / 13,688cm <sup>3</sup>

Carton Dimensions	22 <sup>1</sup> / <sub>4</sub> " x 17 <sup>1</sup> / <sub>2</sub> " x 19"
	56.5cm x 44.5cm x 48.3cm

Shipping Weight	34 lb./15.4kg.
Unit Pack: 1	Cube: 4.28

<b>Optional No. 2180 Support Base</b>	<b>2" / 5cm high</b>
Outside Diameter	12 <sup>3</sup> / <sub>4</sub> " / 32.4cm
Inside Diameter	12 <sup>3</sup> / <sub>8</sub> " / 31.4cm
Elevates Sterilizer Above Surface	1 <sup>3</sup> / <sub>8</sub> " / 3.5cm



## Model 75X (41 qt/39 liter)

**Model 75X-120:** 120 Volt, 50/60 Hz 1650 watts/13.75 amps

**Model 75X-240:** 240 Volt, 50/60 Hz 1650 watts/6.88 amps

<b>Gross Capacity</b>	<b>41 qt / 39 liter</b>
Overall Height	19" / 48.3cm
Bottom Height	14 <sup>1</sup> / <sub>4</sub> " / 36.2cm
Inside Diameter	15 <sup>1</sup> / <sub>4</sub> " / 38.7cm
Unit Weight	45 lbs. / 20.4kg.

### Inner Container No. 4156

Height	10 <sup>1</sup> / <sub>4</sub> " / 26cm
Diameter	14" / 35.6cm
Circumference	44 <sup>1</sup> / <sub>2</sub> " / 113cm
Capacity	27.3 qt / 25.8 liter
Volume	1578in <sup>3</sup> / 25,856cm <sup>3</sup>

Carton Dimensions	24" x 24" x 21"
	61cm x 61cm x 53.3cm

Shipping Weight	51 lb. / 23.1kg.
Unit Pack: 1	Cube: 7

<b>Optional No. 4180 Support Base</b>	<b>3" / 7.6cm high</b>
Outside Diameter	16 <sup>1</sup> / <sub>4</sub> " / 41.3cm
Inside Diameter	15" / 38.1cm
Elevates Sterilizer Above Surface	2 <sup>1</sup> / <sub>4</sub> " / 5.7cm