

FALCON SERIES



SA-300VF

FEATURES:

*** WIDE RANGE OF USAGE :**

"FALCON" autoclaves provide sterilization temperature from 118°C - 134°C. Which are suitable for wrapped / unwrapped materials, utensils, optional(liquids), and etc. Widely used in hospitals, clinics, and laboratories. VERTICAL TYPE SA-300VF are special designed for laboratories, clinics and hospitals.

*** OVER HEAT(LOW WATER) ALARM :**

The over heat will have the buzzer, and stop all the function at the same time. When the water in the chamber is not enough.

*** PRE-HEATING PROGRAM :**

Pre-heating by dry heater, to keep the chamber temperature in balance. To offer the perfect dry function, and the best sterilization function.

SPECIFICATION :

MODEL NO.	SA-300VF
CHAMBER SIZE	Stainless Steel #304 300mmø X 710mm(D)
OVERALL	600(W) X 1090(H) X 450(D)mm
POWER CONSUMPTION	220 - 240V / 2.9KW / 13A
USING TEMPERATURE / PRESSURE	*UNWRAPPED: 118-134°C/2.1kg/cm ² *WRAPPED: 118-134°C/2.1kg/cm ²
DESIGN TEMP.	140°C
SAFETY DEVICES	Pressure control switch, Door safety holder, Over heat protection, Safety Valve, Electronic circuit safety system, Emergency exhaust switch.
STANDARD ACCESSORIES	Stainless Steel #304 Sterilization Basket X 2pc, Heater Cover X 1pc, Clip X 1pc, 2000cc bottle X 1 pc, Exhaust hose X 1pc.
SAFETY APPROVAL	EMC, LVD, by TUV CHAMBER DESIGNED BY ASME/USA STANDARD
AUTO-DRY FUNCTION	YES
CHAMBER CAPACITY	APPROX. 50 LITER